


Highlights of the monotype Z5m implant

Z5m-40-10
> Implant 4mm x10mm



Can be prepped directly in situ:
Abutment and implant shoulder

- individual aesthetics
- Fast and cost-effective restorations – like a natural tooth

- Perfect individual aesthetics due to preppable preparation margin.

Experience from over 60,000 implants placed

- Excellent osseointegration through Z-Systems' proven, patented proprietary SLM® laser surface technology and gentle plasma sterilization
- Self-tapping thread with bone chip reservoir reduces insertion torque

- Conically shaped thread with broadened core diameter in the upper thread area ensures optimal mechanical strength and better primary stability



Z-SYSTEMS
ceramic implants

Z-Systems AG | Werkhofstrasse 5 | CH-4702 Oensingen | support@zsystems.com | www.zsystems.com

Distribución exclusiva: SAT CERAMICS SA. Grecia, Alajuela, COSTA RICA | dental@sat-ceramics.com | www.sat-ceramics.com